

The New
Syro Wave™
Microwave &
Parallel Peptide Synthesizer
- "The Best of Both Worlds"



Biotage Lunch Seminar
5th IPS/47th JPS, Kyoto
Dec 5th 2010



Amit Mehrotra PhD
Marketing Manager – Peptide Chemistry
Product Manager – Synthesis



Biotage

Outline

- Overview of Biotage & MST
- Biotage Peptide Products
 - Peptide Synthesizers
 - ChemMatrix Resins
 - Resolux HPLC Columns
- Microwave Heating
- Syro *Wave*[™]
- Syro *Wave*[™] Application work
- Conclusion

Biotage

With more than 5,000 discovery chemistry systems installed in over 600 facilities worldwide, Biotage automated systems and consumables work together to increase productivity and improve success rates



& PEPTIDES

MultiSynTech

Biotage entered the peptide synthesis business by partnering with MultiSynTech of Germany and provides technical support with fully trained Field Service Engineers and Application Chemists worldwide

- Dr. Udo Treffer, CEO of MultiSynTech, is a peptide chemist, who developed the Syro I and Syro II robotic peptide synthesizers
- Over 150 systems installed
- Over 300 publications in technical journals using Syro systems for peptide synthesis
- Proven reliable in the most demanding applications

Microwave Assisted Peptide Synthesis Using Biotage Instruments

Manual SPPS

M. Erdélyi, A. Gogoll, *Rapid microwave-assisted solid phase peptide synthesis*, *Synthesis*, **2002**, 11, 1592-1596*

M. Brandt, S. Gammeltoft, **K. J. Jensen**, *Microwave heating for solid-phase peptide synthesis: General evaluation and applications to 15-mer phosphopeptides*, *International Journal of Peptide Research and Therapeutics*, **2006**, 12(4), 349-357

Semi-automated SPPS

S. L. Pedersen, K. K. Sørensen, **K. J. Jensen**, *Semi-automated microwave-assisted SPPS: Optimization of protocols and synthesis of difficult sequences*, *Biopolymers (Peptide Science)*, **2010**, 94, 206-212

Fully automated SPPS

L. Malik, A. P. Tofteng, S. L. Pedersen, K. K. Sørensen and **K. J. Jensen**, *Automated 'X-Y' robot for peptide synthesis with microwave heating: Application to difficult peptide sequences and protein domains*, *Journal of Peptide Science*, **2010**, 16, 506-512

* Personal Chemistry

Peptide Market Estimates

- The market for synthetic therapeutic peptides rose from €5.3 billion in 2003 to €8 billion in 2005. It has been estimated that it will reach €11.5 billion in 2013

Pichereau, C. and Allary, C. (2005) Therapeutic peptides under the spotlight. Eur. Biopharm. Rev. 88–91


- CBI (Commonwealth Biotechnologies Inc) believes the market for peptide drugs will achieve a compound annual growth rate of 7.5 per cent and be worth \$13.4bn by 2010 (June 2009)
- The worldwide market for custom peptides is projected to grow at an annual average rate of 11.9%, valuing the market up to \$1bn in 2010

CBI market update, Euroinvestor.co.uk, March 2010

Biotage Peptide Synthesis Systems



ChemMatrix Resins





If your peptide is : long, complex, or hydrophobic: → ChemMatrix resin

- Biotage is now distributing ChemMatrix resins
- Biotage have selected 5 of the most popular linker chemistries for SPPS (Rink, Wang, HMPB, Trityl, PAL)
- ChemMatrix is a patented 100% PEG resin from Matrix Innovation that offers substantial advantages over traditional PS & PEG based resins for SPPS
- Peptides produced with ChemMatrix - higher purity and yields

Benefits of ChemMatrix Resin

- **Exceptional stability** – more stability for the chemistry needed in peptide synthesis
- **No Leaching** – does not add impurity to the end-users work flow
- **Excellent solvent compatibility** – organic or aqueous, water or otherwise
- **Many choices of linker and also pre-loaded options available** – we have access to a wide choice
- **Proven superior performance** – comparison of synthesis of Protease HIV-1 (99 aa) using Wang-ChemMatrix and PS resin shows significant advantages
- **Microwave compatible** – in peptide synthesizers or manual synthesis

Case Study: Protease HIV-1 (99 aa)



After 78 cycles

- **Scale: 0.1 mmole**
- **Conditions: HATU/HOAt/DIEA**
- **Cleavage: Reagent K, 2 h**

Int. J. Pept. Res. Ther. **2007**, 13, 221-227

The Resolux family

- R-P HPLC columns for separation of peptides and proteins
- Analytical, semi-prep & prep columns
- Three product lines, defined on the different pore sizes and each product line is available with two surface chemistries:
 - Resolux 120: A 120 Å media – available with C18 and C8 functional groups
 - Resolux 200: A 200 Å media - available with C18 and C4 functional groups
 - Resolux 300: A 300 Å media - available with C18 and C4 functional groups
- An economical choice

Microwave Heating Gives...

- Faster and more precise heating
- Faster chemical reactions
- Greater yields and better purities
- Novel reactions

Why Does Microwave Heating Speed Up Reactions?


Arrhenius Equation: $k = A e^{-E_a/RT}$

Conventional		Microwave
23°C / 12 h	=	100°C / 5 min
70° C / 16 h	=	150° C / 5 min



Reaction rate roughly doubles for every 10 °C temperature increase

Mechanism of Microwave Heating

Dipolar oscillation



Ionic conduction




Molecular Orientation = Molecular Friction = **Heat**

Advantages of Microwave Heating

- The rate of heating is generally higher than by conventional means
- No direct contact between the energy source and the reaction mixture
 - No temperature gradient through the sample
- The energy transfer is direct to the absorbing reactants
- Allows reactions to occur in a more controlled manner in a decreased time scale


Heating in Synthesis

Indirect heating through vessel



Conventional Heating

Direct heating – transparent vessel



Microwave heating

Syro Wave™ - “The Best of Both Worlds”

This unique Parallel Peptide Synthesizer with Automated Microwave Technology offers “The Best of Both Worlds” for Peptide Chemists



Syro Wave™ - “The Best of Both Worlds”

The Syro Wave is the result of a joint development project announced last year between Biotage and MultiSynTech

It combines the proven performance of the established MultiSynTech robotic synthesizer with Biotage's expertise in microwave technology




Biotage

Syro Wave™ - “The Best of Both Worlds”

- Standalone parallel peptide synthesizers:
 - productivity and cost efficiency
- Standalone microwave peptide synthesizers:
 - difficult or longer peptides
- Inevitably results in:
 - Increased demand on bench space
 - Duplication of computer control systems
 - Multiple operating software platforms to learn and maintain
- Biotage has addressed these shortcomings, with the Syro Wave™, containing both parallel and microwave technologies on the same platform

Biotage Syro Wave™ - Specification



Includes PC, Monitor, Printer
& Syro XP Software

Microwave cavity

Single robotic arm

1 x digital syringe pump

2 Vortex mixers: μ Wave & parallel

24 or 48 position reactor block


Vacuum Pump

Amino Acid Rack

Reagent Bottle Rack

Waste Bottle

Biotage Syro Wave™



Syro Wave™

- The Biotage Syro Wave™ is the **first** valve-free peptide synthesizer with integrated microwave heating
- The system offers unequalled productivity for peptide synthesis with the ability to run 24, 48 or 96 reactions in parallel at RT
- The additional microwave cavity can reduce cycle times and synthesize a long peptide in the shortest possible reaction time
- Unlike other systems that use complex valve modules and compressed air transfer of reagents, the Syro Wave™ system delivers precise volumes of each reagent with a robotic liquid handler. This design minimizes amino acid consumption and therefore reduces synthesis cost
- Proven history of performance and reliability

Syro Wave™ - Key Features & Benefits

- All amino acids and reagents are delivered by a digital syringe pump for the highest accuracy and minimum reagent use
- Dedicated vortex mixer for MW allows for homogeneous heat distribution in the reactor vial
- Disposable reactor vials - avoid wasting time cleaning reactors or replacing blocked frits
- Separate (“off-line”) cleavage & transfer workstation removes any bottlenecks in synthesis



Syro XP Software

The screenshot displays the Syro XP Peptide software interface. The main window features a virtual peptide synthesizer with a 96-well plate and various reagent reservoirs. A menu bar at the top includes File, Activities, Settings, and Help. Below the menu bar are icons for Program, Automatic, and Manual. The main workspace shows a virtual synthesizer with a 96-well plate and several reagent reservoirs labeled DMF, Bottle_5, and 40% Pip/DMF. An 'Activities' window is open, displaying a log of operations. The log includes a 'Stop' button and a message: 'User stopped process at 4:01:02 PM'. The log content is as follows:

```
Activities
[Icons: Fill, Reaction, Empty, Wash, Stop]
Estimated termination on 1/23/2010 at 1:0 ?
Lower 5ml reactor
Cycle1
  MW12_IS_S_10ml_5min75C_exV_start
    Fmoc_Depro_0.1mmol
      FILL - RB_2 [40%Pip/DMF] -
      REACTION - React 02 min; R
      EMPTY - 1; 45 s
      FILL - SL_1 [DMF] -> RV_1
      FILL - RB_2 [40%Pip/DMF] -
      REACTION - React 02 min; R
      EMPTY - 1; 45 s
Execute
Update
User stopped process at 4:01:02 PM
EMPTY
```

Syro XP Software

The Workplace is an active component of the programming desktop



Chemfile - covers all the commands for a reaction step

Workfile - covers all the commands for an entire synthesis cycle

Synthesis file - containing all the cycles for a complete synthesis


Consumables & Accessories




Vial Loading Tool



Consumables & Accessories



5 ml PP reactor vial
with PTFE frit
Disposable



50 ml AA vessel
Re-usable/disposable



Reactor Block Configurations



24 x reactor block

2-ml reactor vials for 5-50 μ mole


or

5 ml reactor vials for 25-150 μ mole

or

10 ml reactor vials for 50-300 μ mole

1 x Syro Wave™ & Syro I, 2 x Syro II



48 x reactor block

2 ml reactor vials for 5-50 μ mole

1 x Syro Wave™ & Syro I, 2 x Syro II




96 x tip reactor (optional)

0.4 ml pipette tips for 1-5 μ mole

1 x Syro Wave™ & Syro I, 3 x Syro II

Cleavage & Transfer Workstation



Parallel cleavage & transfer

3-configurations:

- 24 position reactor block
- 48 position reactor block
- 96 position Tip reactor block

- Pressure mediated transfer of the cleavage solution into the transfer rack
- Simple and reliable manual operation
- This off-line cleavage and transfer process improves throughput

Syro Wave™ - Summary of Advantages

- Microwave and parallel peptide synthesis capability
 - 3 systems in 1
 - Single channel microwave peptide synthesizer
 - Parallel peptide synthesizer RT
 - Single channel peptide synthesizer RT
 - versatile solution to synthesize difficult sequences
 - higher yield for long peptides (and short)
 - faster reactions for reduced synthesis time
 - reduced reagent consumption and lower cost
 - improved reliability

What is a difficult sequence?


So called “**difficult sequences**” are problematic if not impossible to synthesize using standard coupling and deprotection protocols

Difficulties are mainly related to:

- Intra- and/or intermolecular aggregation
- Secondary structure formation
- Steric hindrance of protecting groups which can generate premature termination of the sequence

β -amyloid (1-42) peptide

DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVVIA



β -amyloid (1-42) peptide is the major component of the neurological plaques in Alzheimer's patients

Synthesis is difficult due to reported on-resin aggregation and folding due to hydrophobic C-terminus

β -amyloid (1-42) peptide reagents

- **Resin:** Fmoc-RAM-TG resin 0.24 mmol/g loading
- **Amino Acids:** 770- μ L of 0.5M Fmoc-AA in DMF with HOBt
- **Coupling:** + 880 μ L of HBTU/DMF (0.45 M)
+ 380 μ L of DIPEA/NMP (2.0 M)
- **De-protect:** 2000 μ L of 40% piperidine in DMF
- **Wash:** DMF
- **Notes:** Fmoc-Phe-OH in 0.5M NMP with HOBt (better solubility than in DMF)

β -amyloid (1-42) peptide method

- **Synthesis Scale:** 100 μ mol
- **Deprotection:** 3 min with 40% piperidine in DMF
+10 min with 20% piperidine in DMF
- **Wash:** 4 x 30 sec with DMF
- **Coupling:** 1 x 45 min for RT
or 1 x 5 min @ 75 °C (microwave)
- **Wash:** 3 x 30 sec with DMF

Synthesis of β -amyloid (1-42)

H-DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVVIA-NH₂ (**1**)

Different reaction conditions for the synthesis of **1** using RAM-S-TG resin (0.24 mmol/g):

Entry	Deprotection*	Coupling** Temperature	Coupling time	HPLC purity*** (%)	Overall synthesis time
1	3 + 10 min, RT	RT	1 x 45 min	54	54 h
2	3 + 10 min, RT	75°C	1 x 5 min	73	26 h


*Deprotection: 3 min with 40% pip in DMF followed by 10 min with 20% pip in DMF at RT

**Couplings reagent HBTU / HOBT



***Cleavage: TFA-TES (97.5:2.5) 1ml for 2h. Then add EDT (16 μ l) and TMSBr (13 μ l) for 15 min.

Synthesis of β -amyloid (1-42)

45 min RT, 54% purity



5 min 75 °C, 73% purity



LysM Domain - H-LPERVKVVFPL-NH₂

- The LysM (lysine motif) domain is believed to be involved in the regulation of the interaction between plants and rhizobial bacteria to promote plant growth
- The LysM domain is predicted to consist of two α -helices and a two-stranded anti-parallel β -sheet in a β - α - α - β structure
- The peptide H-LPERVKVVFPL-NH₂ is derived from the C-terminus of the LysM2 domain and is difficult to synthesize as it contains several β -branched and bulky amino acid residues

LysM Domain - H-LPERVKVVFPL-NH₂

H-LPERVKVVFPL-NH₂ (2)

Different reaction conditions for the synthesis of 2 using RAM-S-TG resin (0.24 mmol/g):

Entry	Deprotection*	Coupling** Temperature	Coupling time	HPLC purity*** (%)	Overall synthesis time
1	3 + 10 min,RT	RT	1 x 45 min	52	~14 h
2	3 + 10 min,RT	RT	2 x 45 min	62	~ 24 h
3	3 + 10 min,RT	RT	2 x 120 min	57	~ 52 h
4	3 + 10 min,RT	75°C	1 x 5 min	70	~ 6½ h
5	3 + 10 min,RT	75°C	2 x 5 min	71	~ 9½ h
6	3 + 10 min,RT	75°C	2 x 10 min	77	~ 11½ h

*Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF


**Coupling reagents HBTU/HOBt/DIEA/AA, (3.8:4.0:7.2:4.0)

***Cleavage: TFA-TES-H₂O (95:2:3) 2h.

LysM Domain - H-LPERVKVVFPL-NH₂

2 x 45 min at RT


Purity: 62%



Overall synthesis time: ~24 h

2 x 10 min at 75 °C

Purity: 77%



Overall synthesis time: ~11½ h


1 = product, 2 = des-Leu(N-terminal)

Comparison of Different Reactor Vial Sizes


Synthesis of LysM-Domain C-terminal 11mer: - Coupling: 5 min at 75 °C

- Deprotection: 3 + 10 min at RT


2mL reactor
25 μ mol



5mL reactor
50 μ mol




10mL reactor
100 μ mol




Comparison of Different Reactor Vial Sizes

Synthesis of LysM-Domain C-terminal 11mer: - Coupling: 5 min at 75 °C
- Deprotection: 3 + 10 min at RT


2mL reactor
25 μ mol



5mL reactor
50 μ mol



10mL reactor
100 μ mol



Heating profiles: Blue = temperature
Red = power

Why are we Interested in *N*-Methylated Amino Acids?

- Exist in many biologically-active natural products
- Help obtain information about backbone conformation
- Offer improved lipophilicity, proteolytic stability and bioavailability
- Potentially useful therapeutics

Highly *N*-Methylated Trimer Synthesis

H-MeAla-MeIle-MeGly-NH₂ (**3**)

Different reaction conditions at RT for the synthesis of **3** using RAM-S-TG resin (0.24 mmol/g):

Entry	Coupling Temperature	Coupling time	Coupling reagent*	HPLC purity* (%)	Overall synthesis time
1	RT	20 min	DIC/HOAt	<5	2 h 48 min
2	RT	20 min	HBTU/HOBt	<5	2 h 48 min
3	RT	60 min	DIC/HOAt	<5	4 h 48 min
4	RT	60 min	COMU	<5	4 h 48 min
5	RT	60 min	HATU/HOAt	<5	4 h 48 min
6	RT	2 x 60 min	DIC/HOAt	<5	7 h 51 min
7	RT	24 h	DIC/HOAt	39	~74 h

• Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF.

* Cleavage: TFA-TES-H₂O (95:2:3) 2h.

Highly *N*-Methylated Trimer Synthesis

H-MeAla-MeIle-MeGly-NH₂ (**3**)

Different reaction conditions with microwave heating for the synthesis of 3:

Entry	Coupling Temperature	Coupling time	Coupling reagent*	HPLC purity* (%)	Overall synthesis time
8	75°C	20 min	HBTU/HOBt	8	2 h 48 min
9	75°C	20 min	COMU	59	2 h 48 min
10	75°C	20 min	HATU/HOAt	75	2 h 48 min
11	75°C	20 min	DIC/HOAt	76	2 h 48 min
12	75°C	2 x 10 min	DIC/HOAt	80	3 h 10 min


• Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF.

* Cleavage: TFA-TES-H₂O (95:2:3) 2h.

Highly *N*-Methylated Trimer Synthesis

H-MeAla-MeIle-MeGly-NH₂ (**3**)


24h at RT



Overall synthesis time: ~74 h

Crude Purity
39%

2 x 10 min at 75°C



Overall synthesis time: ~3 h

Crude Purity
80%

ACP-(65-74)

H-⁶⁵Val-⁶⁶Gln-⁶⁷Ala-⁶⁸Ala-⁶⁹Ile-⁷⁰Asp-⁷¹Tyr-⁷²Ile-⁷³Asn-⁷⁴Gly-NH₂ (**4**)

Different reaction conditions using HBTU/HOBt for the synthesis of 4:

Coupling conditions:*

Entry	Residues 65-74	HPLC purity** (%)	Overall synthesis time
1	2 min, RT	85	4 h 53 min
2	2 min, 75°C	90	4 h 53 min


• Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF.

* Coupling reagents HBTU/HOBt/DIEA/AA, (5.0:5.2:10:45.2)

** Cleavage: TFA-TES-H₂O (95:2:3) 2h.


ACP-(65-74)

2 min RT



Overall synthesis time: ~5 h

2 min 75 °C



Overall synthesis time: ~5 h

[⁶⁷Aib-⁶⁸Aib]-ACP (65-74)

H-⁶⁵Val-⁶⁶Gln-⁶⁷Aib-⁶⁸Aib-⁶⁹Ile-⁷⁰Asp-⁷¹Tyr-⁷²Ile-⁷³Asn-⁷⁴Gly-NH₂ (**5**)

Initial screening of coupling reagents and conditions for the synthesis of 5 using RAM-S-TG resin (0.24 mmol/g):

Coupling conditions:

Entry	Residues 66-68	Residues 65, 69-74	Coupling reagent*	HPLC purity** (%)	Overall synthesis time
1	5 min 75°C	2 min 75°C	HATU	0	5 h 43 min
2	5 min 75°C	2 min 75°C	COMU	0	5 h 43 min
3	5 min 75°C	2 min 75°C	DIC	13	5 h 43 min
4	2 x 5 min 75°C	2 min 75°C	DIC	22	5 h 58 min
5	10 min 75°C	2 min 75°C	DIC	20	6 h 19 min
6	2 x 10 min 75°C	2 min 75°C	DIC	30	6 h 49 min

• Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF.

* HOAt was used as additive (except for COMU) and DIEA as base. **3 equiv amino acids.**

** Cleavage: TFA-TES-H₂O (95:2:3) 2h.



Biotage

[⁶⁷Aib-⁶⁸Aib]-ACP (65-74)

H-⁶⁵Val-⁶⁶Gln-⁶⁷Aib-⁶⁸Aib-⁶⁹Ile-⁷⁰Asp-⁷¹Tyr-⁷²Ile-⁷³Asn-⁷⁴Gly-NH₂ (**5**)

Different reaction conditions using DIC/HOAt for the synthesis of 5:

Coupling conditions:*

Entry	Residues 66-68	Residue 65, 69-74	HPLC purity** (%)	Overall synthesis time
7	2 x 10 min, 75°C	2 min 75°C	59	6 h 49 min
8	3 x 10 min, 75°C	2 min 75°C	63	7 h 39 min
9	3 x 90 min, RT	10 min RT	42	20 h 35 min
10	3 x 90 min, RT	45 min RT	37	25 h 50 min
11	3 x 10 min, 75°C	45 min RT	43	13 h 50 min


• Deprotection: 3 min with 40%pip in DMF followed by 10 min with 20%pip in DMF.

* DIC/ HOAt / DIEA. **6 equiv amino acids.**


** Cleavage: TFA-TES-H₂O (95:2:3) 2h.

[⁶⁷Aib-⁶⁸Aib]-ACP (65-74)


3 x 10min 75 °C/ 2min 75 °C



3 x 90 min RT/ 45 min RT



3 x 10min 75 °C/ 45 min RT



1 = Product, 2 = des-Aib

Conclusion: Benefits of Microwave Irradiation in SPPS

- Powerful technique for accelerating the synthesis of peptides and peptidomimetics
- Allows access to the synthesis of difficult peptide sequences
- Improve coupling rates and prevent side reactions
- Use of standard coupling reagents
- Reduction in cost of synthesis
- “Green Chemistry” potential, less harmful solvents

Syro Wave™

Thank you!



Syro Wave™

**We look forward to meeting you on
our Stand**




Biotage