

# Extraction of Mephedrone From Biological Fluids Using EVOLUTE<sup>®</sup> CX

## Introduction

This application note describes the extraction of Mephedrone from biological fluids using EVOLUTE CX.




Figure 1. Structure of Mephedrone

This method demonstrates the extraction of Mephedrone (MCAT) from a range of biological fluids including plasma, urine and whole blood. This extraction is achieved through cation exchange using EVOLUTE CX, a mixed-mode resin-based SPE sorbent. Analytes were quantified using LC-MS/MS and recoveries were all above 90%.

EVOLUTE CX mixed-mode resin-based SPE sorbent extracts a wide range of basic analytes from biological fluids and other aqueous matrices using a generic procedure which minimizes method development time. EVOLUTE CX removes matrix components such as proteins, salts, non-ionizable interferences and phospholipids, delivering cleaner extracts with reproducible recoveries for accurate quantitation.

## Analytes

Mephedrone.

## EVOLUTE CX Configuration

EVOLUTE CX 25 mg 96-well plate, part number 601-0025-P01.

## EVOLUTE CX Procedure

<b>Sample pre-treatment:</b>	Dilute 100 $\mu$ L of sample (plasma, urine, whole blood) 1:3 (v/v) with 50 mM ammonium acetate buffer (pH 6).
<b>Whole blood pre-treatment:</b>	Sonicate 1 mL blood sample for 10 minutes in 50 mM ammonium acetate buffer (pH 6) followed by centrifugation at 5,500 rpm for 10 minutes. Retain supernatant and discard cellular debris (pellet).
<b>Condition:</b>	Methanol (1 mL).
<b>Equilibration:</b>	50 mM ammonium acetate (pH 6) (1 mL).
<b>Sample loading:</b>	Pre-treated sample (400 $\mu$ L).
<b>Interference wash 1:</b>	50 mM ammonium acetate (pH 6) (1 mL).
<b>Interference wash 2:</b>	Methanol (1 mL).
<b>Analyte elution:</b>	5% ammonium hydroxide in methanol (1 mL).
<b>Post Extraction:</b>	Evaporate eluate to dryness and reconstitute in 200 $\mu$ L of water/methanol (80/20 (v/v)) prior to analysis.

## HPLC Conditions

**Instrument:** Waters Acquity UPLC (Waters Assoc., Milford, MA, USA).

**Column:** Acquity UPLC BEH C18 column (1.7 $\mu$ m, 50 x 2.1 mm id) (Waters Assoc., Milford, MA, USA).

**Mobile Phase:** **A:** 0.1% formic acid aq and **B:** 0.1% formic acid/MeOH at a flow rate of 0.5 mL/min.

**Gradient:** Isocratic 60% **A** and 40% **B**.

**Injection Volume:** 5  $\mu$ L.

**Column Temperature:** 35  $^{\circ}$ C.

## Mass Spectrometry Conditions

**Instrument:** Premier XE triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) equipped with an electrospray interface for mass analysis.

**Desolvation Temperature:** 450  $^{\circ}$ C

**Ion Source Temperature:** 150  $^{\circ}$ C

**Collision Gas Pressure:** 3.46 x 10<sup>-3</sup> mbar

Scan Function	Analyte	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
1	Mephedrone (quant)	178.1 > 160.0	35	12
2	Mephedrone (qual)	178.1 > 145.0	35	19

Table 1. MRM transitions for Mephedrone

## Results

All results show recoveries above 90% with RSDs below 10%, Figure 2. shows an example mass chromatogram of Mephedrone (10 ng/mL) from whole blood.




Figure 2. Example mass chromatogram showing Mephedrone at 10 ng/mL from whole blood.

## References

This application note is based on the poster 'Mephedrone: Evaluation of Extraction using Mixed-Mode Cation Exchange SPE with UPLC-MS/MS Analysis.', R Jones et al, presented at TIAFT, Bonn, Germany, Aug 31—Sep 2, 2010.

## Ordering information

Part number	Description	Quantity
601-0025-P01	EVOLUTE CX 25 mg 96-well plate	1

### NORTH AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
press (4) at the auto attendant  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
1-pointsupport@biotage.com

### EUROPE

Main Office: +46 18 56 5900  
Fax: +46 18 59 1922  
Order Tel: +46 18 56 57 10  
Order Fax: +46 18 56 57 05  
order@eu.biotage.com

### Japan

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com

### Distributors

To locate a distributor please visit our  
Web site at [www.biotage.com](http://www.biotage.com)